

NTI-tss PlusTM


Universal Therapy Sets

NDX[®] NATIONAL DENTEX LABS

For more information,
please call 800.678.4140




NTI-tss Plus Universal Therapy Set



NTI-tss Plus Terminator Therapy Set

Simple Steps to a Maximum Protrusive Measurement

- Instruct the patient to protrude their mandible as far as possible
- While in the most protrusive position, use a ruler to measure from the labio-incisal of upper centrals to lingual-incisal of lower incisors



Universal Therapy Set

Dual arch appliances worn to eliminate any tooth contact and accommodate patients with contraindications on both arches. A standard nighttime NTI is worn in conjunction with a flat plane universal slider on the opposing arch.

Terminator Therapy Set

Stainless steel bars are added to the occlusal of the Universal Therapy Set, dramatically increasing the durability. For increased stability and durability, our default is to extend both appliances to the canines.

Indications

- Incisal height difference on both arches
- Patient divots or breaks NTI-tss Plus
- Veneers on both arches
- Missing incisor on both arches

Contraindications

- If concerned about additional vertical opening, Therapy Set may not be recommended.
- Some patients are averse to wearing dual arch appliances.

Patient Care Instructions

- Clean with warm water and anti-bacterial soap after each use
- Do not use alcohol or alcohol-based products for soaking or cleaning
- Never boil the NTI-tss Plus

What to Send

- Full arch VPS impressions, models or digital files
- Maximum protrusive measurement

Insurance Codes

All plans are different so the coverage may vary. In general, the same insurance codes for splints, TMJ appliances, habit appliances all apply to the NTI-tss Plus.

Some useful codes are:

- D9940: Occlusal Guard, by report
- D7880: Occlusal Orthotic Device, by report
- D8210: Removable Habit Appliance

For Medical Submission:

- CPT 21110-52 Interdental Fixation Device/D 7880 Orthotic Appliance
- 2015/16 ICD-10-CM M26.60 Temporomandibular joint disorder, unspecified
- 2015/16 ICD-10-CM M26.69 Other specified disorders of temporomandibular joint

Turnaround Times

Please allow 2 weeks for case turnaround.

NTI TSS
YOUR MIGRAINE GUARD