

## Indications

- Posterior crowns and bridges
- An aesthetic alternative to posterior full-cast crowns and bridges, and PFMs with metal occlusals
- Suitable for implants and posterior crowns and bridges

## Contraindications

- Not desirable for anterior aesthetics

## Benefits

- Made with biocompatible medical-grade zirconia
- Monolithic material makes the restoration nearly chip-proof
- Fabricated using CAD/CAM technology for consistent fit
- Superior strength offering at >1,100 MPa
- Improved gradation of color as compared to other zirconia brands

## Preparation

- Any margin preparation acceptable; slight chamfer preferred
- Conservative preparation similar to a full cast gold
- 1.0-mm of occlusal reduction is ideal.

## Cementing Recommendations

- Cement with a Resin or Resin Modified Glass Ionomer

## Adjusting and Polishing

- Use a fine grit diamond with copious amounts of water and air spray to keep the restoration cool and to avoid micro-fractures.
- *Always polish post-adjustment with porcelain polishing system*

## What to Send

- Upper and lower models, impressions or digital scans
- Shade
- Bite Registration

## Insurance Codes

- D2740- Crown - Porcelain/Ceramic Substrate
- D6245- Pontic Porcelain/Ceramic
- D6740- Abutment Crown Porcelain/Ceramic

**All Verotek restorations have a lifetime guarantee.**