

Designed with hardware on its exterior, dreamTAP[®] offers greater tongue space with no obstructions. The single point of adjustment prevents uneven bilateral adjustment that may create an irregular bite and jaw discomfort. The hardware in dreamTAP[®] is made from nickel-free cobalt chromium.

Features & Benefits

- **Reduces Sleep Apnea:** clinical data shows dreamTAP[®] is the most effective oral appliance available for the treatment of Obstructive Sleep Apnea.¹
- **Successful for severe snoring cases:** AHI (apnea hypopnea index) decreased by 89% in severe Obstructive Sleep Apneatics.²
- **Improves Patient's Sleep:** The quality of their sleep has a huge impact on the overall quality of your patient's life
- Increases lateral and protrusive range of motion for better patient comfort and compliance

Indications

- Anti-snoring
- Obstructive Sleep Apnea
- CPAP Intolerance

Contraindications

- Loose teeth or dentures
- Central Sleep Apnea

What to Send

- Upper and lower VPS impressions or digital files
- Centric bite registration
- Maximum protrusive measurement

Medical Insurance Codes


- G47.33 Diagnostic Code: Obstructive Sleep Apnea
- E0486: Oral Appliance to Treat Sleep Apnea (Medicare Approved)

If you are going to file for insurance to cover the costs of the dreamTAP[®], file for medical insurance, not dental insurance.

These codes have all had successful prior approval. Be sure to have a note of medical necessity signed by the prescriber and a description of the dreamTAP[®] and its function.

¹Thornton Meta-Analysis, 2015. Unpublished

²Hoekema A. et al, Simulated driving in obstructive sleep apnoea-hypopnoea; effects of oral appliances and continuous positive airway pressure (Sleep Breath, 2007), 11:129-138.



Simple Steps to a Maximum Protrusive Measurement

- Instruct the patient to protrude their mandible as far as possible
- While in the most protrusive position, use a ruler to measure from the labio-incisal of upper centrals to lingual-incisal of lower incisors