

To send a file:

- 1 Create an account on www.customer.connectcasecenter.com and follow the prompts for activating the account (if you already have an account, please proceed to step 3).
- 2 After account activation, log into your account and add **your preferred NDX lab** to the list of *Favorite Recipients*.
- 3 Log into your *Primescan Connect* software and complete the digital scan.
- 4 After completing scan acquisition and marking the margin, proceed to the *Connect* step and enter login information.
- 5 Review your order and restoration information, and then click *Enter Order Data*.
- 6 Once the upload is complete, select **your preferred NDX lab** and the desired return date.
- 7 Click *Add Additional Information* and ensure *Anonymize Patient* is not checked unless this feature is wanted for the case. Then, enter gender, age, comments and upload any additional files for the case. If wanting a modelless case, please specify.
- 8 Review information, and then click *Submit Cart*.
- 9 Enter your user name and password once again, and then case submission is complete.



For questions, please contact Kelly Bevington, RDA, EFDA at kelly.bevington@nationaldentex.com or call her at 724.244.9499.