

Carestream® CS Scanner

Offices with a Carestream CS Scanner can send digital impressions directly to NDX nSequence through CS Connect or through direct e-mail.

To send a file:

- 1 Complete the intraoral scan and open the patient's case in the DIS software. This software can operate in stand-alone mode or through Practiceworks.
- 2 Through the file button, in the upper right corner of the patient history window, select *History Open* by clicking the *Book Icon*.
- 3 After highlighting the 3DIO folder, click the *Send To* button located toward the lower right corner of the window. The STL files will now be saved to the desktop.

To submit to NDX nSequence Using CS CONNECT:

- 1 After highlighting the appropriate 3DIO file from step 3, click the *CS Connect* button. The web portal will open to the login screen.
- 2 Log into CS Connect.
- 3 From the drop-down menu, select **NDX nSequence**.
- 4 Attach files by clicking *Browse* and selecting the .dcm file from the patient's folder; then click *Open*.
- 5 Add any pertinent information regarding the case to the comment box.
- 6 Accept the terms and conditions by checking the box and clicking *Submit*.

For questions, please contact your NDX nSequence Sales Specialist.

