

3M™ True Definition Scanner

Doctors with a 3M™ True Definition scanner can send digital impressions directly to NDX Dental Art through 3M Connection Center.

To send a file:

- 1) Log in to your 3M Connection Center Account that was created when your scanner was installed. If assistance is needed, contact 3M Digital Impression Customer Service at 877.722.6528.
- 2) Verify that NDX Dental Art and NDX Dental Art (No Model), are available under the Scanner Laboratory menu.
 - If not, contact 3M Digital Impression Customer Service (877.722.6528) and request adding the profiles to the scanner menu.
- 3) Select the NDX Dental Art profile from the drop down menu that matches your case preference:
 - NDX Dental Art (choose this option when requesting models)
 - NDX Dental Art (No Model)* (choose this option when requesting model-less restorations)
- 4) Complete the intraoral scan and select the NDX Dental Art profile that meets your model preference.
- 5) Complete the True Definition online Rx.
- 6) Upload the scan data to the 3M Communication Center.

* Dentist may request a model-less monolithic restoration only by selecting the NDX Dental Art (No-Model) profile that has been previously added to the scanner menu (see step 2 above).



NDX[®]  Dental Art
YOUR DENTAL LAB PARTNER

3M[™] True Definition Scanner

Should you have questions regarding the products, services and/or processes associated with submitting digital files to NDX Dental Art, please contact us directly at **800.444.3744** or **di-DentalArt@NationalDentex.com**.