

3M™ True Definition Scanner

Doctors with a 3M™ True Definition scanner can send digital impressions directly to NDX Lords through 3M Connection Center.

To send a file:

- 1) Log in to your 3M Connection Center Account that was created when your scanner was installed. If assistance is needed, contact 3M Digital Impression Customer Service at 877.722.6528.
- 2) Verify that NDX Lords and NDX Lords (No Model), are available under the Scanner Laboratory menu.
 - If not, contact 3M Digital Impression Customer Service (877.722.6528) and request adding the profiles to the scanner menu.
- 3) Select the NDX Lords profile from the drop down menu that matches your case preference:
 - NDX Lords (choose this option when requesting models)
 - NDX Lords (No Model)* (choose this option when requesting model-less restorations)
- 4) Complete the intraoral scan and select the NDX Lords profile that meets your model preference.
- 5) Complete the True Definition online Rx.
- 6) Upload the scan data to the 3M Communication Center.

* Dentist may request a model-less monolithic restoration only by selecting the NDX Lords (No-Model) profile that has been previously added to the scanner menu (see step 2 above).



NDX[®] Lords
YOUR DENTAL LAB PARTNER

3M[™] True Definition Scanner

Should you have questions regarding the products, services and/or processes associated with submitting digital files to NDX Lords, please contact us directly at **800.821.0859** or **di-Lords@NationalDentex.com**.