

Midmark Mobile True Definition™ Scanner

Doctors with a Midmark Mobile True Definition™ scanner can send digital impressions directly to NDX Oral Arts Iowa through the Midmark Connection Center.

To send a file:

- 1) Log in to your Midmark Connection Center Account that was created when your scanner was installed. If assistance is needed, please contact Midmark Digital Impression Customer Service at 937.526.8387.
- 2) Verify that NDX Oral Arts Iowa and NDX Oral Arts Iowa (No Model) are available under the Scanner Laboratory menu.
 - If not, please contact Midmark Digital Impression Customer Service (937.526.8387) and request that these profiles be added to the scanner menu.
- 3) Select the NDX Oral Arts Iowa profile from the drop-down menu that matches your case preference:
 - NDX Oral Arts Iowa (choose this option when requesting models)
 - NDX Oral Arts Iowa (No Model*) (choose this option when requesting modeless restorations)
- 4) Complete the intraoral scan and select the NDX Oral Arts Iowa profile that meets your model preference.
- 5) Complete the True Definition online Rx.
- 6) Upload the scan data to the Midmark Communication Center.

* Dentists may request a modeless monolithic restoration by selecting the NDX Oral Arts Iowa (No Model) profile that has been previously added to the scanner menu (see step 2 above).



NDX[®]  Oral Arts Iowa

Midmark Mobile True Definition[™] Scanner

Should you have questions regarding the products, services and/or processes associated with submitting digital files to NDX Oral Arts Iowa, please contact us directly at **800.747.3522** or **di-oralartsia@nationaldentex.com**