

# Carestream CS Scanner

Doctors with a Carestream CS Scanner can send digital impressions directly to NDX Rogers through CS Connect or Direct E-mail.

## To send a file:

- 1) Complete the intraoral scan and open the patient's case in the DIS software. This software can operate in standalone mode or through Practiceworks.
- 2) Through the file button, in the upper right corner of the patient history window, select *history open* by clicking the *book icon*.
- 3) After highlighting the 3DIO folder, click the *Send to* button located toward the lower right corner of the window. The STL files will now be saved to the desktop.

## To submit to NDX Rogers using CSCONNECT:

- a. After highlighting the appropriate 3DIO File from step 3, click the *CS connect* button. The web portal will open to login screen.
- b. Log into CS Connect.
- c. From the drop down menu, select *NDX Rogers*.

Note: If this is the first case being sent to NDX Rogers enter **di-Rogers@NationalDentex.com** in the email field. NDX Rogers will appear in the drop-down menu moving forward.

- d. Attach files by clicking *Browse* and selecting the *.dcm file* from the patient's folder; then click *Open*.
- e. Add any pertinent information regarding the case to the comment box.
- f. Accept the terms and conditions by checking the box and click *Submit*.

A vertical Carestream CS Scanner is shown on the left side of the image. It has a black handle with a green circular button and a white upper section. The text 'CS 3500' is visible on the handle. A dark blue horizontal band is overlaid across the middle of the scanner.

# Carestream CS Scanner

Should you have questions regarding the products, services and/or processes associated with submitting digital files to NDX Rogers, please contact us directly at **800.278.6046** or **di-Rogers@NationalDentex.com**.