

3M™ True Definition Scanner

Doctors with a 3M True Definition scanner can send digital impressions directly to NDX nSequence through the 3M Connection Center.

To send a file:

- 1) Log in to your 3M Connection Center Account that was created when your scanner was installed. If assistance is needed, please contact 3M Digital Impression Customer Service at 877.722.6528.
- 2) Verify that NDX nSequence and NDX nSequence (No Model) are available under the Scanner Laboratory menu.
 - If not, please contact 3M Digital Impression Customer Service (877.722.6528) and request that these profiles be added to the scanner menu.
- 3) Select the NDX nSequence profile from the drop-down menu that matches your case preference:
 - NDX nSequence (choose this option when requesting models)
 - NDX nSequence (No Model*) (choose this option when requesting modeless restorations)
- 4) Complete the intraoral scan and select the NDX nSequence profile that meets your model preference.
- 5) Complete the True Definition online Rx.
- 6) Upload the scan data to the 3M Communication Center.

* Dentists may request a modeless monolithic restoration by selecting the NDX nSequence (No Model) profile that has been previously added to the scanner menu (see step 2 above).



NDX[®]  nSequence

3M[™] True Definition Scanner

Should you have questions regarding the products, services and/or processes associated with submitting digital files to NDX nSequence, please contact us directly at **888.809.2777** or **wetransfer@sequencedental.com**