

3M™ True Definition Scanner

Doctors with a 3M™ True Definition scanner can send digital impressions directly to NDX Green through 3M Connection Center.

To send a file:

- 1) Log in to your 3M Connection Center Account that was created when your scanner was installed. If assistance is needed, contact 3M Digital Impression Customer Service at 877.722.6528.
- 2) Verify that NDX Green and NDX Green (No Model), are available under the Scanner Laboratory menu.
 - If not, contact 3M Digital Impression Customer Service (877.722.6528) and request adding the profiles to the scanner menu.
- 3) Select the NDX Green profile from the drop down menu that matches your case preference:
 - NDX Green (choose this option when requesting models)
 - NDX Green (No Model)* (choose this option when requesting model-less restorations)
- 4) Complete the intraoral scan and select the NDX Green profile that meets your model preference.
- 5) Complete the True Definition online Rx.
- 6) Upload the scan data to the 3M Communication Center.

* Dentist may request a model-less monolithic restoration only by selecting the NDX Green (No-Model) profile that has been previously added to the scanner menu (see step 2 above).



NDX®  Green
YOUR DENTAL LAB PARTNER

3M™ True Definition Scanner

Should you have questions regarding the products, services and/or processes associated with submitting digital files to NDX Green, please contact us directly at **800.247.1365** or digitalscans@greendentalab.com.