

## Indications

- Metal-free or aesthetic RPD required
- When a flexible partial is contraindicated by a flat ridge, a long free end saddle, or few remaining teeth

## Contraindications

- Class II or Class III mobile teeth
- Short clinical crown length or lack of adequate undercut on abutment teeth
- Acrylic allergies
- Limited inter-occlusal space, less than 5mm

## Benefits of ClearFrame

- Metal free aesthetics and increased patient comfort
- Vertical stability with occlusal rests
- Can reline, repair, and polish acrylic saddles chairside
- Less expensive than traditional metal partials
- Hygienic, will not absorb stains and odors
- No metal taste

## Material Properties

ClearFrame<sup>®</sup> is made from state-of-the-art injection-molded nylon material that is strong and non-porous yet easy to adjust and re-polish. The clasps work in the same manner as on a metal framework, engaging undercut on the buccal or labial of the patient's teeth. The partial does not contain Bisphenol A.

## What to Send

- Upper and lower models or impressions with tooth preparations
- Bite registration
- Tooth shade

## Insurance Codes

- D5225 maxillary
- D5226 mandibular